 **SBMS | ULTRA SBMS** **24**
sana bio molecular systems

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**Real-time qPCR instrument
ultraSBMS24**

The text contains the patent data protected by copyright, with all rights reserved. Any part in the document shall never be copied, remade or translated into other languages without the prior written approval of MEDsan® Biotech GmbH.

The analyser is a kind of analyser for special purpose. Please operate in strict accordance with the Operating Instructions to avoid any damage or incorrect results to the analyser.

The user should operate the instrument in strict accordance with the Operation Manual. The protection provided by the instrument may be damaged if not used as specified.

Thank to you for choosing our products. Please make sure to read the operations manual carefully before you use the instrument for the first time!.

Software: SBMSnodus

MEDsan GmbH | Brandshofer Deich 10 | 20539 Hamburg | Germany

After-sales Service: MEDsan® Biotech GmbH

E-mail: devices.service@medsanbiotech.com

Website: <https://www.medsanbiotech.com>

! Statement !

The users should use the instrument strictly in accordance with the method specified in the “Operation Manual”.

If it is not used according to the regulations, the protection provided by the instrument may be damaged.

Contents

CHAPTER 1.	Safety Warning	5
1.1	Use attention	5
1.2	Safety Symbols	8
1.3	Safety precautions	8
1.4	Warning and safety tips	9
1.5	Notice	10
CHAPTER 2.	Introduction	11
2.1	Real-time qPCR instrument	11
2.2	Structural features and module description	11
2.3	Parameters and performance indicators	12
CHAPTER 3.	Instrument installation	12
3.1	Unpack the instrument	12
3.2	Installation and use environment	13
3.3	Power requirements	13
CHAPTER 4.	Instrument use	14
4.1	Boot	14
4.2	Shutdown	14
CHAPTER 5.	Software installation	14
5.1	Windows desktop software installation	14
CHAPTER 6.	Software use	14
6.1	Use of stand-alone operating software	14
6.2	PC operating software use	40
CHAPTER 7.	Maintenance	40
CHAPTER 8.	Troubleshooting	40
8.1	Failure to boot	40
8.2	Software system failure	40
8.3	Test process failure	41
8.4	Abnormal test results	41
8.5	Other faults	41
CHAPTER 9.	After-sales service	41

CHAPTER 1. Safety Warning

1.1 Use attention

In order to use this detection system correctly and effectively, please read the following precautions carefully.

SYSTEM PURPOSE

WARNING

Please refer to «Scope of Application» for the purpose of this detection system. If you need to use this system beyond this scope of application, please consult MEDsan Biotech GmbH. When making clinical judgments based on analysis results, please also consider clinical symptoms or other test results.

OPERATOR

WARNING

This system is limited to inspection professionals or experimenters who have been trained in ultraSBMS24 operation methods to operate and use.

INSTALLATION ENVIRONMENT

WARNING

Please install the system correctly in accordance with the installation environment specified in this manual. Installation and use of this system outside the specified conditions may result in unreliable results and may cause damage to the system.

If you need to change the system status, please contact the customer service department of MEDsan Biotech GmbH or the distributor in your area.

USE ENVIRONMENT

CAUTION:

- (1) Please use the detection system in accordance with the relevant instructions in this manual. Incorrect use may cause incorrect measurement results and may even cause damage to the instrument or personal injury.
- (2) Do not open the upper cover while the detection system is working

ELECTROMAGNETIC COMPATIBILITY

NOTICE

(1) The real-time qPCR instrument (ultraSBMS24 IVD) meets the emission and immunity requirements specified in this part of IEC61326-2-6, see the table below.

Table 1:

Electromagnetic emission	
Launch test	Compliance
CISPR 11 Conducted emission	1 group A
CISPR 11 Radiation emission	
IEC61000-3-2 Harmonic emission	Not applicable
IEC61000-3-3 Voltage fluctuation/flicker emission	Not applicable

Table 2:

Electromagnetic immunity			
Immunity test item	Basic standard	Test value	Meet performance
Electrostatic discharge (ESD)	IEC61000-4-2	Contact discharge: ±2kV, ±4kV Air discharge: ±2kV, ±4kV, ±8kV	B
Radio frequency electromagnetic field	IEC61000-4-3	3V/m, 80MHz ~ 2.0GHz, 80% AM	A
Burst	IEC61000-4-4	Power cord: ±1kV (5/50ns, 5kHz)	B
surge	IEC61000-4-5	Line to ground: ±2kV Wire to wire: ±1kV	B
RF conduction	IEC61000-4-6	Power cord: 3V/m, 150kHz ~ 80MHz, 80%AM	A
Power frequency magnetic field	IEC61000-4-8	3A/m, 50Hz/60Hz	A
Voltage dips, interruptions	IEC61000-4-11	1 cycle 0%; 5/6 cycle 40%; 70% for 25/30 cycle; 250/300 cycle 5%	B C C C

Performance judgment:

- A. During the test, the performance is normal within the specification limits.
 - B. During the test, the function or performance is temporarily reduced or lost, but it can recover by itself.
 - C. During the test, the function or performance is temporarily reduced or lost, but operator intervention or system reset is required
-
- (2) The user is responsible for ensuring the electromagnetic compatibility environment of the equipment so that the equipment can work normally.
 - (3) It is recommended to evaluate the electromagnetic environment before using the equipment.

WARNING

- (1) The **ultraSBMS24** IVD real-time qPCR instrument is designed and tested according to the Class A equipment in CISPR 11. In a home environment, this device may cause radio interference, and protective measures are required.
- (2) It is forbidden to use this equipment near strong radiation sources (such as unshielded radio frequency sources), otherwise it may interfere with the normal operation of the equipment.

1.2 Safety Symbols

 **CAUTION** and  **WARNING**

“Caution”, “Warning” and “Notice” are all mentioned in this manual.

Types of safety signs: The following signs appear in this manual or instrument accessories:

 **CAUTION** — If you do not follow the prescribed procedures or instructions, it may cause damage to the equipment.

 **WARNING**—If you do not follow the prescribed procedures or instructions, it may cause physical injury.

1.3 Safety precautions

Prevent electric shock

To prevent electric shock, please observe the following precautions.

WARNING

When the main power of the detection system is turned on, unauthorized maintenance personnel must not open the fixed casing of the instrument. Spilled reagents and samples inside the instrument may cause the instrument to malfunction and cause electric shock.

Don't place samples and reagents above the instrument. If spilling occurs, turn off the power immediately and contact the service department of MEDsan GmbH | Brandshofer Deich 10 | 20539 Hamburg | Germany.

Protection against biological hazards

In order to effectively protect against biochemical hazards, please observe the following precautions.

Biological infection risk:

If you use infectious samples, incorrect use of samples may result in infection. Do not touch the sample mixtures and waste liquids directly with your hands. Please be sure to wear gloves and overalls during operation to prevent infection and wear protective glasses if necessary.

If the sample accidentally touches the skin, please handle it immediately in accordance with laboratory work standards and consult a doctor.

Prevent injury from static electricity

To prevent injury caused by static-charged parts, please observe the following precautions.

Beware of static electricity:

There may be static electricity on the surface, please do not touch it during use.

Prevent injury caused by hot parts

To prevent burns, please observe the following precautions.

Warning:

Here show fever or high temperature, do not touch during use.

1.4 Warning and safety tips

Please read the following warnings carefully before use and must be strictly followed. Warning: Please read the following precautions carefully before using this instrument

- Avoid blood, reagents, staples and other metal parts from entering the inside of the instrument.
- If an abnormality occurs, cut off the power immediately and unplug the power plug from the power socket. At this time, an inspection application should be submitted to the distributor or our company's agent immediately.
- The operator should not touch the electronic circuit in the instrument, especially contact with wet hands is more likely to get an electric shock.
- When repairing and inspecting the instrument, you must wear rubber gloves and use the specified tools and parts. After the operation, please wash your hands with disinfectant. Otherwise, the part of the skin in contact with blood may be infected or get an electric shock or be burned.
- Be very careful when handling specimens. Be sure to wear rubber gloves, otherwise it may cause infection. If the specimen enters the eye or wound, rinse it immediately with plenty of water and accept a doctor's examination.
- If the equipment is not used in accordance with the method specified by this instrument, the protection provided by the instrument may be destroyed.

Use and disposal of reagents (If used)

- Avoid contact of reagents with skin and clothing during operation.
- If the reagent gets into the eyes accidentally, rinse it off with plenty of water immediately and accept a doctor's examination.
- If you drink the reagent by mistake, you should immediately seek help from a doctor and drink plenty of water to spit out the reagent.
- If the reagent gets on your hands or skin, immediately rinse it off with clean water.
- Waste products such as used reagent strips and other equipment consumables should be treated as medical waste or infectious waste. If it is contaminated by blood, it may be infected by pathogens.

Power supply voltage, connection and grounding

- Ensure that the network power supply and grounding environment of the instrument is good and stable.
- Do not insert the power plug into a power outlet outside the allowable voltage range.
- When installing the instrument, be sure to use the three-core power cable that comes with the machine and ensure that it is well grounded and placed in a position where it is easy to operate and disconnect the power supply. Otherwise, it may cause fire or electric shock.
- When connecting peripheral equipment, be sure to cut off the power first. Otherwise, it may cause electric shock or malfunction.

1.5 Notice

Please pay attention to the following precautions during use:

- (1) When a failure occurs, you can refer to the «Troubleshooting» section for exclusion. If the emergence of the problem cannot be solved, please contact our service department.
- (2) It is forbidden to use this equipment near strong radiation sources (such as unshielded RF sources), otherwise it may interfere with the normal operation of the equipment.
- (3) When using this real-time qPCR instrument, in order to get the best PCR detection results, it is recommended to use paraffin oil.

CHAPTER 2. Introduction

2.1 Real-time qPCR instrument

Product name: Real-time qPCR instrument Model specification: ultraSBMS24 IVD

Scope of application: Medical institutions, testing institutions, food companies, breeding companies, schools and other research and teaching markets.

Intended use: This product is based on fluorescent polymerase chain reaction, which can analyze nucleic acid samples (RNA/DNA) derived from the human body in medical and clinical applications.

Transportation conditions and methods: The product should be placed upwards during transportation to prevent heavy pressure, rain and direct sunlight, and the maximum number of stacked layers cannot exceed 2 layers. Transport according to the storage temperature in 3.2.

The real-time qPCR instrument is a detection instrument completely independently developed by our company, which can complete the gene amplification of 24 holes and 4 channels.

Support online update of control software functions, support USB and WIFI control modes.

2.2 Structural features and module description

The equipment frame is shown in Figure 2.2.1:


Figure 2.2.1

Module description:

Hot lid: Make sure the temperature above the test tube is higher than the temperature of the test tube itself.

Heating module: Control its temperature to raise and lower temperature.

Temperature sensor: Collect the temperature of the thermal lid or heating module to control the output.

Refrigeration module: For air cooling, cooling the heating module. Optical acquisition module: Data acquisition during PCR reaction. Light source: Fluorescence excitation.

Main control module: Control all peripherals for temperature rise and fall, optical acquisition, temperature acquisition, etc.; At the same time, it communicates with the human-computer interaction module.

Human-computer interaction module: Control screen, communicate with main control module and computer.

Computer: Communicate with human-computer interaction module via USB/WIFI.

2.3 Parameters and performance indicators

Parameters	Performance indicators
size	235 * 212 *255mm
Maximum heating rate	≥5.0 °C/s
Maximum cooling rate	≥3.0 °C/s
Module temperature control accuracy	≤0.2 °C/s
Module temperature uniformity	±0.25 °C/s
Temperature control range	35 - 99°C
Detection throughput	24

CHAPTER 3. Instrument installation

3.1 Unpack the instrument

3.1.1 Unpack step

Unpack the instrument and remove the pearl cotton and other materials used for transportation. Save the packing box and packing materials so that you can repack the instrument in the future. The specific steps are as follows:

- (1) Remove the equipment
- (2) Take out the accessories and check whether the items are complete according to the accessories list. If there are any missing, please notify the service apartment or the dealer immediately.

Note: The product packing accessories are subject to the packing list. If any parts are missing or the accessories do not match the packing list, please contact the seller.

3.1.2 Transport method

When the machine is moved out of the packing box, pay attention to keeping the machine stable.

Note: Please save the packing box for packing and transportation before long-distance transportation. The instrument must be placed on a horizontal operating table, not on an inclined surface.

The installation process of the instrument is handled by professionals authorized by our company. To ensure the normal operation of the instrument, please choose an indoor workplace that meets the following requirements to place the instrument:

- The ground or table top is flat;
- No direct sunlight;
- No large amount of dust;
- No strong electromagnetic radiation;
- Easy to operate the location of the disconnect device;
- Well-ventilated location;
- Avoid humidity and high temperature. Avoid violent vibration and shock.

3.2 Installation and use environment

- Normal working environment of the instrument is 15 °C ~ 30 °C , relative humidity 20~70% , no condensation.
- The storage environment of the instrument is: the packaged UltraSBMS24 IVD instrument should be stored at -20°C ~ 55°C, relative humidity≤85%, no corrosive gas and a well-ventilated environment.
- After the instrument is installed, try to avoid frequent movement, and avoid vibration if it needs to be moved.
- Must be installed and moved by authorized professionals.
- When installing, make sure that the distance between the side of the instrument and other objects for quality inspection is not less than 20cm to facilitate ventilation and heat dissipation.

3.3 Power requirements

- Power supply voltage: 100V ~ 240V
- Frequency: 50Hz / 60Hz
- Power: 500 W

CHAPTER 4. Instrument use

4.1 Boot

In order to use the real-time qPCR instrument safely, correctly and stably, you must complete the boot process in the prescribed order:

- (1) Install the antenna to the fixed position of the antenna behind the instrument;
- (2) Plug in the power cord;
- (3) Turn on the power switch on the back of the device.

4.2 Shutdown

- (1) Turn off the power switch on the back of the device.

CHAPTER 5. Software installation

5.1 Windows desktop software installation

- (1) Go to www.medsanbiotech.com → support Download the software and example files.
- (2) A 64 -bit system is required to install the software.
- (3) If you have any questions about the installation, please contact customer service staff to consult the installation process.

CHAPTER 6. Software use

6.1 Use of stand-alone operating software

6.1.1 Main interface Profile

After entering the operating software, the software enters the opening animation, and then the main interface is displayed as shown in Figure 6.1.1.1. At this time, the top cover of the device is in a locked state, click the lock icon in the lower right corner to unlock and lock the top cover:

Figure 6.1.1.1

The interface function is divided into six operation units: [New Detect], [Login], [Setup], [Load Template], [Records], [Lock Lid].

Click the label 1 in Figure 6.1.1.1 to edit the serial number of the device, as shown in Figure 6.1.1.2: The serial number can be used to identify your devices in case you operate more devices. This serial number is not equal the device technical serial number.


Figure 6.1.1.2

After entering the device number, click [Save] to edit the device number. The device serial number is used to identify the corresponding device that needs to be connected when the PC is connected.

Click the Bluetooth icon  marked 2 in Figure 6.1.1.1 to connect to the Bluetooth printer, and a list of Bluetooth devices will pop up in this interface, as shown in Figure 6.1.1.3:


Figure 6.1.1.3

Click the Bluetooth printer to be connected for pairing connection, as shown in Figure 6.1.1.3.1:


Figure 6.1.1.3.1

After clicking [Ensure], the pairing code will pop up (according to the pairing code in the manual of the paired Bluetooth printer, generally 0000 or 1234), as shown in Figure 6.1.1.3.2:


Figure 6.1.1.3.2

After the pairing and connection are completed, the connected state of the printer will be displayed. After the binding, the operation does not need to be repeated next time when restarting, as shown in Figure 6.1.1.3.3:


Figure 6.1.1.3.3

Click the network icon  labeled 3 in Figure 6.1.1.1, the wireless network list pops up in this interface, as shown in Figure 6.1.1.4:


Figure 6.1.1.4

Select the network you want to connect to connect. After the network connection is successful, click the red cross button to exit the interface.

Click the mark 4 in Figure 6.1.1.1 to perform the system time setting function. After the time is selected, click the [Ensure] button to complete the system time setting, as shown in Figure 6.1.1.5:


Figure 6.1.1.5

6.1.2 New Detect

After opening the lid of the instrument and putting in the reagents, click the main interface function [New Detect] to enter the function interface of [Create Detect] setting, as shown in Figure 6.1.2.1:


Figure 6.1.2.1

The user name in this operation interface is not editable. After logging in, the user name of the account will be displayed. Enter the detect name and select the detect type.

There are two types of detect: QPCR and Standard. Select [QPCR] or [Standard], and click [Next]. Enter the function interface of operating condition setting as shown in Figure 6.1.2.2 and Figure 6.1.2.3. Figure 6.1.2.2 is for the setting of QPCR operating conditions, and Figure 6.1.2.3 is for the setting of Standard operating conditions:


Figure 6.1.2.2

QPCR condition setting:

As shown in figure 6.1.2.2, in the condition input box, you can slide left and right to view or edit the conditions, click the icon marked **1** [+] in the figure to set the temperature conditions, click the mark **2** [-] in the figure to delete the selection To set the temperature conditions, click the mark **3** in the figure to enter the number of cycles that need to be run, click mark **4** in the figure to enter the temperature, click mark **5** in the figure to enter the time, click mark **6** in the figure to choose whether to collect, click mark **7** in the figure to enter the reagent volume, click mark **8** in the figure [**Save template**]. click mark **9** to add 1-3 PCR cycles, click mark **10** to delete 1 PCR cycle.

Standard setting: as shown in Figure 6.1.2.3, enter the temperature and time to complete the input of the conditions, click [**Next**] to enter the sample setting function, as shown in Figure 6.1.2.4:

Figure 6.1.2.3

Click **[Next]** to enter the sample setting function, as shown in Figure 6.1.2.4:

Figure 6.1.2.4

Select the hole to be tested, click **[Editor]** to edit the sample information of the hole, this function can edit multiple holes at the same time, as shown in Figure 6.1.2.5:

Figure 6.1.2.5

Enter the sample name, test item, and barcode in this interface (the three parameters can be omitted), select the fluorescent dye (required), select the Task flag information for each channel (Unknow, negative, standard), click **[Ensure]** to complete the sample setting, and click **[Run]** to start the test.

When prompted like in Figure 6.1.2.6, please confirm that the lid has been closed, and click [Ensure] to start the test:


Figure 6.1.2.6

After the device is started, it will jump to the running status interface, as shown in Figure 6.1.2.7:

Figure 6.1.2.7

In this interface, the progress animation, lapsed time, number of cycles run and estimated remaining time are displayed. Click the **[Stop]** button to stop the test, and slide to the left to switch the operating condition setting interface, including the display animation of the running status, which will display the current running stage of the instrument in real time, as shown in Figure 6.1.2.8:


Figure 6.1.2.8

Continue to slide to the left to display the collected data curve, as shown in Figure 6.1.2.9:

Figure 6.1.2.9

Click the button marked 1 in the figure to display all the curves and select the hole position to display the data curve on the hole position; after the test is completed, the test completed window will pop up on the interface and display the function button of [Determine], such as Figure 6.1.2.10.


Figure 6.1.2.10

Click the [Determine] function button to enter the result analysis interface, as shown in Figure 6.1.2.11:

Figure 6.1.2.11

In Figure 6.1.2.11, the result interpretation uses a list to display the result data, label **1** is the channel to be analyzed, label **2** is the modified threshold. Ct value will be calculated after the threshold is entered. Click [**Save**] to store the analyzed results and click the red cross in the right corner to return to the Run curve interface.

Details on template settings, channel selection and result interpretation as well as result determination can be found at:

6.1.3 Log in, registration, password modification and other functions

Click the [**Login**] button in the main interface to jump to the login interface, this function can be used to log in, as shown in Figure 6.1.3.1:

Figure 6.1.3.1

Enter the account password to log in (**Note:** the account password needs to be registered by the salesperson in the background, after logging in, the result data can be uploaded to the cloud service, which can be viewed on the PC side)

Click [**Forgot Password**] to modify the account password, as shown in Figure 6.1.3.2. Please enter the user name, new password, administrator email account, administrator email password, and then click [Send] to send the verification code to the administrator email, and enter it in the interface, and click [**CHANGE PASSWORD**] to retrieve the password.

Figure 6.1.3.2

Change password

Click the [Change Password] button to enter the password modification interface. After entering the account number, the original password, and the new password, click the [CHANGE PASSWORD] button to modify the password.

Registration

Click the [REGISTER] button to enter the registration interface, select the [Administrator] item through the selection, first register the account as an administrator, then fill in the corresponding account information, complete the verification of the mailbox, and complete the registration;

After the administrator is registered, the user can select the [User] item through the selection, and enter the account and password to register as a common user.

SBMSnodus

Enter Username	⊗
Enter Password	🗑
Administrator	▼
Email account	@gmail.com
Email password	
Enter Verification Code	SEND

REGISTER

Figure 6.1.3.3

6.1.4 Setup

Click the [Setup] button in the main interface to jump to the operation interface of function settings, in which export record, result print, language setting, software update, and device information query functions can be performed, as shown in Figure 6.1.4.1:

Operator	Detect name	Time	Detect type
user	2022_01_14_07_28_19	2022/01/14/ 07:28:45	QPCR
user	ultraSBMS SARS-CoV-2-I	2022/01/14/ 06:47:35	QPCR
user	2022_01_14_06_44_38	2022/01/14/ 06:45:09	QPCR
user	2022_01_14_10_08_46	2022/01/14/ 10:09:21	QPCR

Figure 6.1.4.1

Click **[Export Records]** to enter the operation interface of record management, as shown in Figure 6.1.4.1, click the **[Management]** button to perform record management, as shown in Figure 6.1.4.2:

Figure 6.1.4.2

Select the record to be deleted, exported to USB in Excel, sent to email in excel, printed in A4 paper or printed by thermal print. There are two export options, you can select to export EXCEL. Whichever way you choose, you need to insert a USB flash drive. After the USB flash drive is inserted, the interface will display the USB connection logo, as shown in Figure 6.1.4.3:

Figure 6.1.4.3

Click **[excel]** to export the original data file in xls format. By selecting email or usdisk, it can be exported to the corresponding mailbox or inserted USB flash drive as shown in Figure 6.1.4.4:


Figure 6.1.4.4

By entering the sender's email address, sender's password and the recipient's mailbox, the data can be sent by email.


Figure 6.1.4.4.1

The exported data is stored in the “ultraSBMS24 IVD” folder in the root directory of the USB flash drive.

After inserting the USB flash drive, the USB flash drive needs to be ejected safely, otherwise the USB flash drive may be damaged or the exported data may be incomplete.

Click USB connection icon  on the interface, then the logo prompts whether to eject the USB stick. After clicking [Ensure], it prompts that the eject is successful, and the USB flash drive can be ejected safely, as shown in Figure 6.1.4.5 and Figure 6.1.4.6:


Figure 6.1.4.5


Figure 6.1.4.6

Click the [A4 Print] button turn to the printing and edit page of PDF. The edit interface will pop up as shown in Figure 6.1.4.7

Data Check
Inspection Report
X

Name

Surname

Sex

ID Card / ID No.

Date of Birth

Date of Sampling

Phone

Email

ZIP Code City

Sample Type

Street Name Hous Nr

Hole	Project	Ct	Result
1	SARS-COV-2-II(COVID-19)	-	NEGATIVE
2	SARS-COV-2-II(COVID-19)	-	NEGATIVE
3	SARS-COV-2-II(COVID-19)	-	NEGATIVE
4	SARS-COV-2-II(COVID-19)	-	NEGATIVE
5	SARS-COV-2-II(COVID-19)	-	NEGATIVE
6	SARS-COV-2-II(COVID-19)	-	NEGATIVE

Medical director/operator

Figure 6.1.4.7

You can print in the PDF editing interface. Click the [Print] button to jump to the system print interface, as shown in Figure 6.1.4.8. Manually select a network printer under the same local area network, select the page to be printed and set related printing information, and click the print icon on the upper right, as shown in Figure 6.1.4.9, to print the file.


Figure 6.1.4.8


Figure 6.1.4.9

Click [**Project result**], and enter the surface as shown in figure 6.1.4.10.

Figure 6.1.4.10

Click **[Manage]** button to perform record management, as shown in Figure 6.1.4.11. The result of the project can be printed by A4 paper or thermal print paper:

Figure 6.1.4.11

Select the results to be printed. This function needs to be connected to a Bluetooth printer. For connection details, please refer to the operation in Figure 6.1.1.3. Click **[Print]** to pop up a confirmation window, click **[Ensure]** to print the result.

The export and print PDF functions are the same as the **[Run Record]** function.

You can find more about the report function at:

Click **[Reagent Project]** to add reagent items and parameter setting functions, as shown in Figure 6.1.4.13


Figure 6.1.4.13

In this function interface, click **[Add]** to add reagent items; as shown in Figure 6.1.4.14, in this interface, pass Scan the reagent QR code on the reagent kit to add data, or click **[Manually]** to manually enter the parameters to enter the reagent parameters. After the data is entered, click **[Save]** to complete the addition of reagent items.

Figure 6.1.4.14

Click the displayed item, You can view the detailed parameters of this project, as shown in Figure 6.1.4.15. In this interface, click [Editor] to edit the parameters of the reagent project. After editing, click [Save] to save the data parameters.

Figure 6.1.4.15

Click [Manage], then you can delete the reagent item, as shown in Figure 6.1.4.16. Select the reagent item that needs to be deleted, and click [Delete] to delete the selected reagent item


Figure 6.1.4.16

Click [**Language Setting**] to set the German and English language settings, as shown in Figure 6.1.4.17. Select the language to be displayed, and the software will restart itself

Figure 6.1.4.17

Click [**Software Updates**] to update the software online, as shown in Figure 6.1.4.18. After clicking update, the software will automatically update and restart:


Figure 6.1.4.18

Click [**About Instrument**] to enter the display interface of device information, as shown in Figure 6.1.4.19:

Figure 6.1.4.19

6.1.5 Load the template

Click [**Load the Template**] in the main interface to enter the template management interface, as shown in Figure 6.1.5.1:


Figure 6.1.5.1

Click the list option to directly enter the function interface of the detect settings. For details, please refer to 6.1.2 for detect settings. Click the **[Manage]** button to organize the templates as shown in Figure 6.1.5.2. Note that the built-in templates in the system cannot be deleted.

Figure 6.1.5.2

6.1.6 Records

Click **[Records]** in the main interface to enter all test records under the current user and non-login, as shown in Figure 6.1.6.1

Figure 6.1.6.1

Click the list option to directly enter the function interface of the experiment setting. For details, please refer to 6.1.2 for experiment setting.

Connect to PC control:

This function needs to be connected to the PC through Type-C or connected to the PC under the same WIFI; please see Figure 6.1.1.4 for the connection steps for network connection.

6.2 PC operating software use

You can find availability and more infos here:

CHAPTER 7. Maintenance

The user needs to clean the surface of the equipment and the turntable regularly to keep the equipment tidy. In the event of dumping or leakage of reagents or samples, carry out anti-pollution treatment according to the relevant operations in 1.4. No cleaning agents or disinfectants that are incompatible with the equipment parts or the materials contained in the equipment have not been found. If you have any questions about the compatibility, please contact the after-sales department.

Do not use the device when it is being repaired or there is a fault that needs to be dealt with.

If it is not used for a long time, please pack the equipment, and place it in a dry place.

CHAPTER 8. Troubleshooting

8.1 Failure to boot

- (1) The main power switch is not turned on. →→→ Turn on the main power switch.
- (2) The power cord is not connected properly. →→→ Please reconnect the power cord.
- (3) Others. →→→ Please contact the after-sales department.

8.2 Software system failure

- (1) System operation failure. →→→ Please contact the after-sales department.
- (2) The test analysis software cannot be started. →→→ Please contact the after-sales department.

8.3 Test process failure

8.3.1 The display is not working

- (1) White screen, blurred screen or no response to touch. →→→ Please contact the after-sales department.

8.3.2 The test process stops suddenly

- (1) It may be that the power supply of the tester is interrupted. →→→ Turn on the tester switch again and test again.
- (2) A communication failure occurs. →→→ Turn on the tester switch again and test again.
- (3) The problem still exists. →→→ Please contact the after-sales department.

8.4 Abnormal test results

- (1) The fluorescence curve is abnormal. →→→ The brightness of the light is abnormal, check the tightness of the instrument; check whether the experiment parameter settings are normal; and check whether the sample and reagent are contaminated.
- (2) Re-test, if the problem still exists. →→→ please contact the after-sales department.

8.5 Other faults

When other failures occur, please contact the after-sales service department in time!

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