MEDsan® Transfer Set

INTENDED USE

This product is used for the collection, pretreatment, and release of nucleic acids of SARS-CoV-2 samples. It is intended to be used in combination with MEDsan SARS-CoV-2 RT PCR kits for the detection of SARS-CoV-2 RNA. It is suitable for the qualitative RNA detection from SARS-CoV-2 samples in nasal, oropharyngeal (throat), oral/ saliva, combined nasal/oropharyngeal, anterior nasal, or mid-turbinate nasal swabs based on the type of the swab from individuals suspected of COVID-19 by their healthcare provider or for screening of individuals without symptoms or other reasons to suspect COVID-19 infection.

Detection in combination with a validated RT PCR kit is achieved via rapid nucleic acid extraction from respiratory tract samples followed by multiplex real-time RT-PCR targeting SARS-CoV-2 in real-time PCR instruments.

PRINCIPLE OF PROCEDURE

The swabs can be taken from different locations and can be a combination of two extraction locations. Nasal, oropharyngeal (throat), oral / saliva, combined nasal/oropharyngeal, anterior nasal, or mid-turbinate nasal swabs can be used for the sample collection based on the type of the swab.

COMPONENTS

There are 2 versions of this product:

1 x 100 Tests (T02001)

2 x 100 Tests (T02002)

Please check the contents on the sticker of the package

T02001	Transfer & Extraction Tubes	2x50 tubes (2 ml)
	Specimen Collection Swabs	1x100
T02002	Transfer & Extraction Tubes	2x50 tubes (2 ml)
	Specimen Collection Swabs	2x100

Note: The different batches are not recommended to be mixed.

STORAGE RECOMMENDATIONS AND EXPIRATION DATE

- Validity: 12 months.
- Storage Temperature (Preservation solution tubes): 2°C ~ 30°C
- The production date and expiration date are shown in the outer package.
- Do not use any part of the package if the expiration date has passed.
- A damaged package or a damaged component of the package should no longer be used.

SAMPLE REQUIREMENTS

1. Specimen Type:

Nasal swabs, oropharyngeal (throat) swabs, oral/ saliva swabs, combined nasal / oropharyngeal swabs, anterior nasal swabs, mid-turbinate nasal swabs.

2. Collection of Clinical Specimens:


- Make sure you are wearing your personal protective equipment and gear before you start collecting samples (gloves, eye protection, face mask, protective clothing etc.).
- Open the sealed tubes bag and prepare one tube for each swab.
- Determine the type of sampling you want to proceed to.
- Remove the swab from the package, and keep the swab tip clean.
- Take a sample from the desired sampling location.
- While ensuring the swab does not come into contact with anything, carefully remove the cap of the tube.

- Immediately place the swab inside the specimen tube. Do not remove the liquid from the specimen tube. Break the swab at the breakpoint and close the cap tightly.
- Leave the swab in the buffer solution for at least 5 minutes. After 5 minutes, the sample solution is ready for real-time PCR using the compatible RT PCR kits and thermocyclers. If the RT PCR is not performed immediately, store the samples at 2-8°C.
- Place the specimen container in a biohazard bag or store it in the fridge at 2-8°C.
- Change your gloves before taking a sample from another person.

For more information refer to the CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19.


Nasal specimen collection

- Ask the patient to take off the mask and to blow their nose to clear nasal passage of excessive mucus.
- Tilt the patient's head back 70 degrees.
- Carefully insert a swab into the nostril until resistance is met.
- Gently rotate the swab several times against the nasal wall to absorb secretions. Repeat in the other nostril using the same swab.
- Place the swab into the extraction tube containing the buffer solution.
- Ask the patient to wear a mask.

Oral saliva specimen collection

- Ask the patient to take off the mask.
- Ask the patient to rinse the mouth with clean water before sampling.
- Ask the patient to cough up secretions deep in the throat and accumulate in the oral cavity above the tongue.
- Insert the specimen collection swab into the mouth and place it over the tongue.
- Rotate the swab to absorb the sample for 10-15 seconds and remove it.
- Place the swab into the extraction tube containing the buffer solution.
- Ask the patient to wear a mask again.

Oropharyngeal specimen collection

- Ask the patient to take off the mask.
- Tilt the patient's head back 70 degrees
- Take the single use sterile specimen swab. Insert the swab into the mouth to the posterior pharyngeal and tonsillar areas and gently rotate it. Avoid touching the tongue and teeth.
- Place the swab into the prepared extraction tube.
- Ask the patient to wear a mask again.

3. Storage and Shipment of Clinical Specimens:

Specimens must be packaged, shipped, and transported according to the current edition of the International Air Transport Association (IATA) Dangerous Goods Regulation.

The sample solution in the test tube can be stored at 4°C for up to 14 days.

LIMITATIONS OF THE TEST

Unreasonable collection, transportation, and storage of specimens, improper operation, and improper laboratory environment may lead to inaccurate results.




If the stored specimen is used after a storage period longer than 24 hours, false-negative results cannot be excluded.

This test has been authorized only for the detection of SARS-CoV-2, not for any other pathogens.

PRECAUTIONS

- For *In vitro* diagnostic use only.
- For single use only. Do not reuse.
- Do not perform the test in a room with strong airflow and in an environment that is too hot, too humid, or too dry.
- Handle all specimens as if infectious by using safe laboratory procedures.
- Do not mix components from different lots.
- Make sure to correctly label your specimens.
- Wearing protective gear is obligatory for the whole process of specimen collection.
- The liquid is non-toxic and harmless. If it splashes on the skin, wipe the buffer and wash the skin with water.

INDEX OF SYMBOLS

Symbol	Meaning
	<i>In-vitro</i> diagnostic medical device
	Use by date
	Temperature Limit
	Keep away from sunlight
	Keep dry
	Do not use if package is damaged
	Lot Number/Batch Code
	Manufacturer
	Sufficient for <n> tests
	Consult instructions for Use
	European Conformity

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MEDsan® GmbH