


Guangzhou Four E's Scientific Co., Ltd

5th floor, Building C, No. 2, Ruitai Road,
Huangpu District, Guangzhou, 510535, Guangdong, China
Tel: 86-20-3803 2775; Fax: 86-20-3803 2771

Single-Channel Adjustable Volume Pipette



Features:

- Ergonomic design of low plunger pressure avoids repetitive strain injuries
- Volume range of 0.1ul-10ml
- Large 4-position volume display and high display accuracy up to 0.002ul
- Durable high-tech synthetic material provides excellent resistance to high temperature, corrosion and weather
- Reasonable structure and exquisite technology make flexible assembling and disassembling.
- Tip ejector makes easy and quick tip change
- High-elastic tip cone provides good seal and easy removal of tips

- Easy to calibrate and maintain with tools supplied
- Different volume range is identified by different color top
- Calibrated in accordance with ISO8655 and each pipette supplied with individual test certificate
- The lower section can be autoclaved

Specifications:

Single-channel Adjustable Volume Pipettes				
Model	Volume Range	Increment	Tip Type	Color
LH0301001	0.1-2.5µl	0.002ul	10ul	Brown
LH0301002	0.5-10µl	0.01ul	10ul	Orange
LH0301003	2-20µl	0.02ul	200ul	Purple
LH0301004	5-50µl	0.05ul	200ul	Green
LH0301005	10-100µl	0.1ul	200ul	Dark Blue
LH0301006	20-200µl	0.2ul	200ul	Blue
LH0301007	100-1000µl	1ul	1000ml	Red
LH0301008	1000-5000µl	5ul	5ml	Yellow
LH0301009	1-10ml	10ul	10ml	Grey